

# Predicting HIV Viral Load


Hyojin Lim and Tom Lindman

---

# HIV Viral Load Testing

---

- > **Without a cure, antiretroviral therapy (ART) is central to treating HIV<sup>1</sup>**
- > **Viral load testing essential to ART**
- > **Viral load testing essential to combating the HIV epidemic**



# Viral Load Testing: Challenges

- > Quantitative polymerase chain reaction (qPCR) the “gold standard” for viral load testing<sup>2</sup>
- > qPCR requires:
  - Expensive lab equipment
  - Highly trained staff
  - Time
- > Implementation challenges in low- and middle-income countries<sup>3</sup>


2. Zhao et. Al (2019)

3. Greig et. Al (2020)



# Viral Load Testing: Potential Solutions

- > Isothermal DNA amplification offers potential qPCR alternative
- > Reaction produces florescent spots on paper slip
  - Spots correlate with viral load
  - Can process tests using phone application
- > Challenge: prediction at high spot concentration
- > Project goal: increase prediction range, accuracy



# User and Use Cases

---

## > Users: Main Project team

- Cole Martin and collaborators (primary users)
- Use cases:
  - > Update phone app to predict viral load using our model
  - > Update our model using new data, features, and ML algorithms

## > Users: Clinicians testing phone app


- Clinicians working with patients on ART
- Use case:
  - > Monitor patient viral load using our model
  - > Interact with model exclusively via phone app



# Prediction with ML models

---

- > ML models are effective in analyzing data
- > Lightweight ML models can be run on mobile devices
- > Key points:
  - Feature Selection: Use minimal, relevant data
  - Model Selection: Opt for lightweight & accurate models.
  - Result Analysis: Evaluate to enhance accuracy.



# Feature selection

---

## > Key features: 11 features in total

- max\_fluorescence\_change:  
Maximum change in fluorescence
- max\_area\_change:  
Maximum area change
- max\_spot\_change:  
Maximum spot change

Critical features for  
extending the model's  
dynamic range

- percent\_area:  
Percentage of area covered
- n\_spots:  
Number of spots identified

Less critical



# Model selection

## > Models tested

- **Linear Regression: Simple, interpretable**
- **Random Forest: Handles non-linearity, reduces overfitting**
- **XGBoost: Fast, high-performance gradient boosting.**



	Linear Regression	Random Forest	XGBoost
<b>Average Train MSE</b>	0.0090	0.0028	0.0000
<b>Average Test MSE</b>	0.0127	0.0152	0.0138





# Result analysis and key takeaways

## > Result analysis



## Future work

---

- > **Collect more data**
- > **Conduct sensitivity analysis**
- > **Fully integrate Android App**

